



Active wheelchair

Model 1.360 / 1.370

Operating manual



MEYRA®
We move people.

Contents

Introduction	5
Important information	6
List of models	6
Indications	7
Acceptance	7
Specifications	8
Use	8
Adjustment	9
Reinstallation	9
Life span	9
Overview	10
Model 1.360	10
Model 1.370	11
Brake	12
Locking the pressure brakes	12
Releasing the pressure brakes	13
Leg supports	14
Calf belt	14
Adjusting the calf belt	14
Removing/attaching the calf belt	14
Lower leg support	15
Footplates	15
Footboard	16
Leg support upper part	17
Turning the leg supports to the side	17
Swivelling in the leg supports	17
Removing the leg supports	18
Attaching the leg supports	18

Clothes guard	19
Back support	20
Back support belt	20
Fitting the back belt	20
Push handles	21
Height adjustable push handles	21
Height adjustment of the push handles	21
Removing the push handles	22
Inserting the push handles	22
Seat	22
Seat belt	22
Utensils bag	22
Wheels	23
Tyre damage on pneumatic tyres	23
Drive wheels	23
Quick release axle	23
Quick release axle with tetra locking device Code 697	24
Support castors	25
Swinging the support castors	25
Folding/Unfolding	26
Folding the wheelchair	26
Carrying the wheelchair	26
Unfolding the wheelchair	26
Transport	27
Transport security	27
Transport of people inside a motor vehicle	27
Maintenance	28
Maintenance	28
Information for the specialist dealer	28
Maintenance schedule	29

Technical data	32
Model 1.360	33
Model 1.370	36
Meaning of the labels on the wheelchair	39
Meaning of the symbols on the type plate	40
Inspection certificate	41
Warranty / Guarantee	42
Warranty / Guarantee section	43
Inspection certificate for transfer	43

INTRODUCTION


We thank you for the confidence you have placed in our company by choosing a wheelchair from this series.

The model of your selection, fulfils the wish for mobility and more independence by way of a new styling of the proven MEYRA technology.

With all equipment and their accessories the wheelchair offers die respective adaptation to your disability.

Additional information about our products can be found on our website under:

< www.meyra.com >.

 Contact your specialist dealer when required.

Like any other vehicle, a wheelchair is a technical aid. It is subject to explanations, requires regular care and can cause danger when used improperly. The correct handling must therefore be learned. This manual is to help you get accustomed to the handling of the wheelchair as well as to prevent accidents.

 **Note:**

Please note that the illustrated equipment variants can deviate from your model.

We have therefore also listed chapters with options that might not be applicable for your individual wheelchair.

Important information

Attention:

- ! Read and observe the following documentation belonging to the wheelchair before first operation:
- this operating manual,
- the safety and general handling instructions < *Mechanical and muscle powered wheelchairs* >.

🗣️ Note:

Children and juveniles should read the documentation belonging to the wheelchair together with their parents respectively a supervisor or an accompanying person before first use.

Users with visual impairments can find the PDF-files of above mentioned documents on our website under:

< www.meyra.com >.

- 🗣️ Contact your specialist dealer when required.

Alternatively users with visual impairments can have the documentation read out by a helper.

LIST OF MODELS

This operating manual applies to the following models:

Model 1.360

Model 1.370

INDICATIONS

If the following indications occur we recommend the application of this mobility product:

- ☞ Walking disability resp. extremely limited walking ability as part of the basic need to move around in your own home.
- ☞ The need to be able to leave home for a short walk in fresh air or in order to reach the places, commonly in the perimeter of the home, required to fulfil basic needs.
- ☞ The ability to use the adaptive wheelchair with own personal strength must be given.

Provision with an adaptive wheelchairs can be an option, when the adaptation and adjustment possibilities of standard or light-weight wheelchairs are not sufficient.

ACCEPTANCE

All products are checked for faults in the factory and packed in special boxes.

☞ **Note:**

However, we request that you check the vehicle for possible transport damage immediately on receipt – preferably in the presence of the carrier.

☞ **Note:**

The packaging of the wheelchair should be stored for a further transport that might become necessary.

SPECIFICATIONS

The wheelchair was developed for adults and adolescents.

The wheelchair solely serves to transport one person in the seat and not as a hauling aid, transporter or similar.

USE

The wheelchair is applicable on level, firm surfaces and can be used as follows:

- for indoors (e.g. apartment, day care),
- outdoors (e.g. paved paths in parks),
- as a companion on tours (e.g. in a bus or train).

The wheelchair offers manifold adjustment possibilities to individual vital statistics.

The optimised crossbrace construction, apart from the small folding dimension, offers high stability.

The wheelchair should be adapted to your needs by a specialist dealer before the first use. The adaptation will take into account the driving experience, the physical limits of the user and the main place of use of the wheelchair.

Attention:

- ! Always have adaptation and adjustment work carried out by a specialist dealer.

ADJUSTMENT

The specialist workshop will hand out the wheelchair to you under consideration of all relevant safety instructions, ready for operation and adjusted to your needs.

🔧 Note:

- We recommend a regular control if the wheelchair adjustment in order to ensure a long-term optimal provision even with changing illness/handicap patterns of the user.
- 🔧 We recommend regular medical exams in order to ensure safety for active participation in traffic.
- 🔧 Retrospective adjustments should be carried out solely by the specialist dealer!

REINSTALLMENT

The wheelchair is suited for reinstallation. With the building block system the wheelchair can be fit to accommodate different handicaps body sizes.

LIFE SPAN

We expect an average life span of about 4 years for this product, as far as the product is applied for its designated purpose and all maintenance and service guidelines.

The life span of your product depends upon the frequency of use, the application environment and care.


The implementation of spare parts can prolong the life span of the product. As a rule spare parts are available up to 5 years after production is discontinued.

- 🔧 The indicated lifespan does not constitute additional guarantee.

OVERVIEW

Model 1.360

The overview shows the most important components of the wheelchair.




Pos. Description

- | | |
|---------------------------------|---------------------------------------|
| (1) Back support | (8) Locking knob / Quick release axle |
| (2) Clothes guard / arm support | (9) Handrim |
| (3) Seat cushion | |
| (4) Footboard | |
| (5) Steering wheel | |
| (6) Driving wheel | |
| (7) Pressure Brake | |

OVERVIEW

Model 1.370

The overview shows the most important components of the wheelchair.



Pos. Description

- | | |
|---------------------------------|---------------------------------------|
| (1) Back support | (8) Locking knob / Quick release axle |
| (2) Clothes guard / arm support | (9) Handrims |
| (3) Seat cushion | |
| (4) Leg support | |
| (5) Steering wheel | |
| (6) Driving wheel | |
| (7) Pressure Brake | |


BRAKE

By locking the pressure brakes [1] with the brake lever (2), the wheelchair is secured against rolling away unintentionally (parking brake).

 **Note:**

Therefore observe the maintenance instructions as well as safety and general handling instructions < *Mechanical and muscle powered wheelchairs* > chapters < *General safety information* > and < *Brakes* >.

Attention:


 Arrange an immediate repair of the brakes by your specialist workshop if the braking performance reduces.

Locking the pressure brakes

To secure the wheelchair against any unintentional rolling, press both brake levers down (3) resp. outward (4) as far as possible.

 **Note:**



It should not be possible to push the wheelchair forward when both brakes are locked.



Releasing the pressure brakes

Swivel both brake levers up (5) resp. inward (6) as far as possible.

- ⚠ Do not grab into the area of the cross-brace [7] when releasing the pressure brake (6). – Danger of jamming the fingers!



LEG SUPPORTS

Attention:

- ! Before any actions on the leg supports the wheelchair is to be secured against unintentional rolling motions. – View chapter *Brake* on page 12

Calf belt

The calf belt (1) prevents the feet from sliding off the back of the footboard.

Adjusting the calf belt

The calf belt is guided through the front frame tubes, view arrows [2], or around special attachment pins [3] and adjusted in length with a velcro strap.

Removing/attaching the calf belt

Removing and attaching the calf belt is achieved with the velcro fastener.

Attention:

- ! Never drive without calf belt (except when scutteling)!
- ☞ < Scuttling > is slow movement of the wheelchair with the feet. For this the footboard is folded up.



Lower leg support

For entry or exiting as well as „scutteling“ (the slow propulsion of the wheelchair with the feet) the footboard [1] / the footplates [2] can be folded upward.

- ☞ Check the locking points!
 - Remove both feet from the footplates.
 - Remove lower calf belt, if present.
 - ☞ Therefore observe chapter *Calf belt* on page 14.

☞ Note:

Before starting to drive the footboard / the footplates need to be folded back down [4] and the calf belt attached.

Footplates

The footplates can be folded outward and up [2] resp. inward and down [4].



Footboard

The footboard [1] can be folded up to one side.

 **Note:**


Fold the footboard up before swivelling away and removing the upper leg support.



Folding up the footboard

In order to fold up the footboard lift the loose end of the footboard as far as possible [2].

Folding the footboard down

In order to fold down the footboard, lower the loose end of the footboard as far as possible down onto the footboard bracket [3].

 Ensure that the footboard has safely locked into place.



Leg support upper part

The upper leg support with an inserted lower leg support is termed leg support.

Turning the leg supports to the side


For easy transfer out of/into the wheelchair as well as driving closer to a closet, bed or bathtub the leg supports can be swivelled away toward the in-/outside (1).

Note:

Before removing the leg supports remove the calf belt and fold the foot-plates up.

-  View chapter *Calf belt* on page 14 and chapter *Lower leg support* on page 15

Attention:

-  Leg supports turned to the side are released automatically and can easily come off. Note this when handling (e.g. transport).


Pull or press the respective locking lever (2) backward and swivel the corresponding leg support outward.

Swivelling in the leg supports

For inward swivelling, let the leg supports swivel forward until the lock audibly engages [3].

Note:

After audibly swivelling the leg supports inward check the respective locking device.

-  Afterwards observe chapter *Lower leg support* on page 15.



Removing the leg supports


For easy transfer into and out of the wheelchair as well as a reduced wheelchair length (important for transport) the leg supports can be removed [1].

Note:

Before removing the leg supports loosen or remove the calf belt on one side.

 View chapter *Calf belt* on page 14

First swivel the leg support sideways [2] and then remove them toward the top.


 Therefore observe chapter *Turning the leg supports to the side* on page 17.

Attaching the leg supports

With the leg support in a swivelled-aside position, hang it in and then swivel it to the front until the locking device audibly latches.

Note:

After attachment swivel the leg supports inward.

 Therefore observe chapter *Swivelling in the leg supports* on page 17.

Press the leg supports, swivelled to the side, parallel to the front frame tube and lower it into place [2]. – In doing so the holding pin must slide into the frame tube.



CLOTHES GUARD

The clothes guard [1] protects the clothes from the open spokes of the drive wheels and from spray water.

Attention:

- ! No not grab between the clothes guard and tyre. – Danger of squashing!
- Do not lift or carry the wheelchair by holding on to the clothes guard / arm supports.
- When the wheelchair is being pushed by an accompanying person the user is to place his hands onto the arm support pads (2) or in his lap and not at the sides between body and clothes guard. – Danger of squashing the fingers!



BACK SUPPORT

Back support belt


Fitting the back belt

The tensioning of the back support [1] is adjustable.

- Pull off the flap of the back support and fold it to the front [2].
- Open the Velcro fastener of the belt (3) that you wish to adjust and close it again after adjustment.

Attention:

- ! The overlapping of the loop and hook strap has to be at least 10 cm!
- Fold the flap of the back support back down and fasten it with the velcro fastener [4].



Push handles

The height adjustable push handles [1] are attached to the back support with a clamping jig each (2) and can be adjusted to the requirements of the accompanying person.

Height adjustable push handles

The push handles are continuously height adjustable and secured against pulling out.



Height adjustment of the push handles

- First hold onto the push handle that is to be adjusted with one hand, then loosen the corresponding clamping lever (2) with the other hand.
- Then lift/slide the push handle to the desired height and afterwards retighten the clamping lever (2).

Special features of the clamping lever

The clamping lever (2) can be turned into an operation position that is comfortable for you.


- For this pull the lever outward, until the teething is released.
- After turning the lever let the teething engage again.



Removing the push handles

Loosen the clamping lever (3) and pull the according push handle up as far as it will go.

Press the respective spring button (4) inward and pull the push handle out of the clamping device (5).



Inserting the push handles

Press down the spring button (4) and insert the push handle from the top into the respective clamping device (5).

Slide the push handle through the clamping device (5) and retighten the clamping lever in the desired position (3).

SEAT



Seat belt

Utensils bag

A utensils bag (1) with velcro fastener is located underneath the seat belt.

In order to open the utensils bag (1) pull the leather latch down (2).

To close the utensils bag press the front edge back up.



WHEELS

Tyre damage on pneumatic tyres


- For repairing tyre damage we recommend the use of a foam cartridge that is available in speciality shops. – Afterwards look up a specialist workshop as soon as possible.



Drive wheels

The drive wheels are on a full quick release axle (1).

- No person may be seated in the wheelchair during assembly or removal. The wheelchair must stand on a level and firm surface. Before starting the disassembly work, support the frame to prevent the wheelchair from tipping over and secure it to prevent an unwanted movement or tipping over.



Note:

After each assembly the locking device is to be checked with sideward pulling/pressing on the drive wheels!

- If the drive wheels has too much sideward lag or the quick release axle does not engage, contact your specialist dealer immediately for repair.

Attention:

- After inserting the drive wheel the locking knob (1) must stick a couple of millimetres out of the axle nut.

Quick release axle

The driving wheels can be removed and re-assembled without any tools.

- First press the locking knob (1) of the quick release axle in the center of the hub.
- Afterwards remove or attach the drive wheel [2].

Quick release axle with tetra locking device Code 697




The driving wheels can be removed and re-assembled without any tools.

Removing the drive wheel:

- First turn the activation wheel (1) of the quick release axle as far as possible about 90° in the direction of the rim [2].
- ☞ The marker (3) faces in the direction of the rim when the quick release axle is released [2].
- Pull off the drive wheel [4].

Attaching the drive wheel:

- If necessary previously turn the activation wheel (1) as far as possible in the direction of the rim [2].
- Slide the axle of the drive wheel as far as possible into the axle receptacle.
- Afterwards turn the activation wheel (1) of the quick release axle by 90° in the direction of the center of the wheel hub (5).
- ☞ The quick release axle is locked when the marker (3) is located centred to the wheel hub (5).



SUPPORT CASTORS

The support castors [1] serve to raise the tilting stability and can be swivelled inward underneath the seat (2).

- ☞ Depending on the version the support castors can be attached singly or as a pair.

Swinging the support castors

Press the support castors out of the locking then swivel them inward until the lock engages automatically [2].



FOLDING/UNFOLDING

Folding the wheelchair

Your wheelchair can be folded in a few moments without tools being required [1].

- Remove the seat cushion, if applicable.
- Remove the calf belt if necessary.
- Remove the leg supports or footboard, resp. fold up both footplates, view chapter *Leg supports* on page 14
- Bend the back belt to the rear.
- Pull the seat belt upward in the centre of the front and back [1].

☞ Ensure that the folding fixation (3), if existent, engages correctly.

Carrying the wheelchair

Your wheelchair can be carried without difficulty when folded.

- From the front push one lower arm under the upward folded seat belt.
- Place other hand underneath the rear seat fold for support.
- Lift the wheelchair to horizontal position.


Unfolding the wheelchair

For unfolding the wheelchair lift it slightly on one side and loosen the folding fixation (3), if existent. At the side resting on the floor, press the seat tube down to the end stop [2].

☞ Note:

You may need to press down both seat tubes with the hands in order to do this.

☞ Observe that the seat profile lies securely on the support surface.



TRANSPORT

Transport security

The wheelchair is only to be secured through the four securing points [1] and (2).



- ☞ The anchor positions are marked with the symbol (3).
- ☞ The procedure for securing the wheelchair can be read in the document Safety and general handling instructions < *Mechanical and muscle powered wheelchairs* > chapter < *Transport in motor vehicles or conveyors* >.

Transport of people inside a motor vehicle

To determine if your wheelchair is approved as a seat for transport inside a motor vehicle, please look at the type plate of your wheelchair.

Attention:

- ! Subsequently, firmly attached components of non-original parts on your wheelchair are not approved for person transportation within a motor vehicle.



MAINTENANCE

An incorrect or neglected cleaning and maintenance results in a limitation of the product liability.

Maintenance

The following maintenance Instruction gives you a guide for carrying out the maintenance work.

- ✎ They do not give information about the actual extent of work required on the vehicle.

Information for the specialist dealer

The maintenance and service manual for this product can be found on our website under:

< www.meyra.com >.

You can for example find the following information in these manuals:

1. Adjustments that can be carried out with tools.
2. Step by step explanations to important repairs.
3. Information on model specific amendments.
4. A checklist for the annual inspection.

The functional tests necessary for the inspection are listed in the check list.

They are a guide for the performance of the inspection work.

✎ Note:

It does not outline the actual scope of the necessary work which can only be ascertained by an inspection of the vehicle.

After the successful completion of an annual inspection the inspection certificate should be recorded in the operating manual.

Maintenance schedule

WHEN	WHAT	Remark
Before starting out	<p>Test brakes for faultless operation</p> <p>Activate brake lever to the limit. The locked wheels should not be able to turn under operating conditions. If they can still turn, the brakes must be repaired by an authorised specialist workshop.</p>	Carry out test yourself or with a helper.
Before starting out	<p>Check pressure brake for wear</p> <p>Move brake lever to the side</p>	<p>Carry out tests yourself or have a helper do it.</p> <p>If you notice any increasing slackness on the brake lever take the wheelchair to your specialist workshop immediately for repairs.</p> <p>– Danger of accident!</p>
Before starting out (when applicable)	Check air pressure of the tyres	<p>Do it yourself or with the aid of a helper.</p> <p>For this use an air pressure gauge or, if not available, „Thumb-pressure-method“ or similar (view < <i>Brakes</i> > in document Safety- and general handling instructions < <i>Mechanical and muscle powered wheelchairs</i> >.</p>
Before starting out Especially before driving in the dark	<p>Check the lighting</p> <p>Check the lighting equipment and reflectors for flawless functioning.</p>	Carry out test yourself or with a helper.

WHEN	WHAT	Remark
Every 2 months (depending on distance covered)	Check frame tubes for damage	Carry out test yourself or with a helper. If deformations or cracks occur in the welding seams, contact a specialist workshop immediately for repairs. – Danger of accidents!
Every 2 months (depending on distance covered)	Lubricate the following components with a few drops of oil – Bearing of crossbrace. – Moving parts of the locking mechanism. – Brake lever bearings.	Do it yourself or with the aid of a helper. Components must be free from used oil residues before lubrication. Please ensure that excess oil does not contaminate the environment (e.g. your clothing).
Every 2 months (depending on distance covered)	Check all screw connections for secure fit	Carry out test yourself or with a helper. The following is to be checked with great care: – Attachment of the back and seat belt, – Attachment of the footboard.
Every 2 months (depending on distance covered)	Check tyre profile Minimum tread = > 1 mm	Carry out a visual check yourself or with a helper. If the tyre profile is worn down or if the tyre is damaged, consult a specialist workshop for repairs.
Every 6 months (depending on frequency of use)	Check: – Cleanness. – General condition.	View < <i>Maintenance</i> > in document Safety- and general handling instructions < <i>Mechanical and muscle powered wheelchairs</i> >.

WHEN	WHAT	Remark
Manufacturer recommendation: Every 12 months (depending on frequency of use)	Maintenance jobs – Vehicle	To be carried out by the specialist dealer.

TECHNICAL DATA

All model dependent specifications within the < *Technical data* > refer to the standard version of the respective models.

The overall length depends on the position and size of the drive wheels.

If not noted otherwise the dimensions are determined with drive wheels of \varnothing 610 mm (24").

The widths were determined with a handrim distance of 15 mm.

Dimensional tolerance ± 1.5 cm, $\pm 2^\circ$.

Short form of wheelchair dimensions:

SH = Seat height

SW = Seat width

SD = Seat depth

BH = Backrest height

Calculation of the max. user weight:

Attention:

- ! The permissible total weight is calculated from the empty weight of the wheelchair and the maximum user- (person-) weight.

Additional weight due to subsequent additions or luggage reduce the maximum permissible passenger weight.

Example:

A driver wishes to take luggage with a weight of 5 kg. Thus, the maximum user weight is reduced by 5 kg.

Tyre filling pressure:

Maximum tyre pressure is printed on the tyres on each side.

Full tyre pressure – steering wheel

Standard:

2.5 - 3.5 bar = 36 - 50 psi

Full tyre pressure – drive wheel

Standard:

3.0 - 4.0 bar = 44 - 58 psi

Ultra-light running tyres:

6 bar = 87 psi

High pressure:

8 bar = 116 psi

Model 1.360

All data within the following table relates to the standard version of the stated model.

Dimensional tolerance ± 1.5 cm, $\pm 2^\circ$.

Model: **1.360**

Type plate: at the crossbrace tube

Life span: 4 years

Dimensions:.....min. / max. / manufacturer setting

Overall length (frame Code 346):

Short frame:840 / 910 / --- mm

Medium frame:880 / 950 / --- mm

Long frame:920 / 990 / --- mm

Overall length (frame Code 347):

Short frame:810 / 880 / --- mm

Medium frame:850 / 920 / --- mm

Long frame:890 / 960 / --- mm

Overall width:

w wheel camber 0.5° :510 / 670 / --- mm

w wheel camber 3° :560 / 720 / --- mm

Overall height:650 / 1000 / --- mm

Back strap height adjustable by +2.5 cm each:

Standard/adjustable back support strap:25 / 50 / --- cm

Seat width:340 / 500 / --- mm

Seat depth (short frame 380 - 420 / medium 440 - 460):380 / 460 / --- mm

Seat depth (long frame):480 / 500 / --- mm

Seat height, without cushion:

Seat surface height at front edge:460 / 520 / --- mm

Height of armrests:

from seat surface:15 / 29 / --- cm

Back support to front edge of arm support:27 / 27 / --- cm

Seat cushion thickness:3 cm

Push handle height Code 502:

continuously adjustable: max. 15 cm

Dimensions:.....min. / max. / manufacturer setting

Back support angle:

bent at manufacturer (firm):..... 84° / 96° / 96°

Seat tilt:..... 0° / 15° / ---

Leg support angle:

Code 346..... 77°

Code 347..... 85°

Foot support to seat, without seat cushion (lower shank length):

with foot board:..... 39 / 48 / --- cm

Wheels

Steering wheel:

ø 100 x 30 mm:.....solid rubber

ø 125 x 30 mm:.....solid rubber

ø 142 x 37 mm:.....solid rubber

Drive wheel (pneumatic tyres):

ø 596 mm (24 x 1"):.....8.0 bar

ø 616 mm (24 x 1 3/8"):.....4.0 bar

ø 596 mm (24 x 1"):.....puncture safe

ø 610 mm (24 x 1 3/8"):.....puncture safe

Handrim diameter:..... 53.5 cm

Wheel camber:..... 0.5° / 3° / ---

Axle horizontal position:

Active:..... 80 / 100 / 80 mm

Middle:..... 50 / 70 / 50 mm

Passive:..... 30 / 50 / 30 mm

Transport dimensions:

Folding length:..... 810 / 910 / --- mm

Length without drive wheels:..... 610 mm

(Support castors are removed or swivelled underneath the seat)

Folding width:

w wheel camber 0.5°:..... 300 mm

w wheel camber 3°:..... 350 mm

Folding height:..... 650 / 1000 / --- mm

Permitted inclination/slopes

max. obstacle height (depending on the setting of the leg support height):..... 0 to 100 mm
min. turning radius diameter:..... 1100 mm
Max. permissible rising gradient:4.5° (8 %)
Max. permissible falling gradient:4.5° (8 %)
Max. permissible transverse gradient:4.5° (8 %)
static tilting safety in all directions:.....6° (10 %)

Climatic data:

Ambient temperature:.....-25 °C to +50 °C
Storage temperature:-40 °C to +65 °C

Weights

Permissible total weight:..... max. 129 kg
max. user weight (including additional load):..... max. 120 kg
Max. additional loading: 10 kg
Empty weight:min. 9 / max. 15 kg
Heaviest single component:
Drive wheel:2.2 kg
Transport weight:.....min. 6.5 kg
(without arm supports, cushion, drive wheels)

Model 1.370

All data within the following table relates to the standard version of the stated model.

Dimensional tolerance ± 1.5 cm, $\pm 2^\circ$.

Model: **1.370**

Type plate: at the crossbrace tube

Life span: 4 years

Dimensions:.....min. / max. / manufacturer setting

Overall length (with leg support):

Short frame:..... 900 / 980 / ---- mm

Medium frame:..... 920 / 1020 / ---- mm

Long frame:..... 960 / 1060 / ---- mm

Extra long frame:..... 1000 / 1100 / ---- mm

Overall length (w/o leg support):

Short frame:..... 720 / 800 / ---- mm

Medium frame:..... 740 / 840 / ---- mm

Long frame:..... 780 / 880 / ---- mm

Extra long frame:..... 820 / 920 / ---- mm

Overall width:

w wheel camber 0.5° :..... 490 / 690 / --- mm

w wheel camber 3° :..... 540 / 740 / --- mm

Overall height: 650 / 1000 / --- mm

Back strap height adjustable by +2.5 cm each:

Standard/adjustable back support strap:..... 25 / 50 / --- cm

Seat width: 320 / 520 / --- mm

Seat depth (short frame): 360 / 400 / --- mm

Seat depth (medium frame):..... 420 / 440 / --- mm

Seat depth (long frame):..... 460 / 480 / --- mm

Seat depth (extra long frame):..... 500 / 520 / --- mm

Seat height, without cushion:

Seat surface height at front edge:..... 460 / 540 / --- mm

Height of armrests:

from seat surface:..... 15 / 29 / --- cm

Back support to front edge of arm support:..... 27 / 27 / --- cm

Seat cushion thickness: 3 cm

Push handle height Code 502:

continuously adjustable: max. 15 cm

Back support angle:

bent at manufacturer (firm): 84° / 96° / 96°

Seat tilt: 0° / 15° / ---

Leg support angle:

Code 93 77°

Foot support to seat, without seat cushion (lower shank length):

with single footplates: 37 / 51 / --- cm

Wheels

Steering wheel:

ø 100 x 30 mm: solid rubber

ø 125 x 30 mm: solid rubber

ø 142 x 37 mm: solid rubber

ø 180 x 40 mm: pneumatic tyres

Drive wheel (pneumatic tyres):

ø 596 mm (24 x 1"): 8.0 bar

ø 616 mm (24 x 1 3/8"): 4.0 bar

ø 646 mm (26 x 1"): 8.0 bar

ø 596 mm (24 x 1"): puncture safe

ø 610 mm (24 x 1 3/8"): puncture safe

Handrim diameter (24"/26"): 53.5/ 58.5 cm

Wheel camber: 0.5° / 3° / ---

Axle horizontal position:

Active: 80 / 100 / 80 mm

Middle: 50 / 70 / 50 mm

Passive: 30 / 50 / 30 mm

Transport dimensions:

Folding length: 880 / 1100 / --- mm

Length without drive wheels, without leg supports: 530 / 850 / --- mm
(Support castors are removed or swivelled underneath the seat)

Folding width:

w wheel camber 0.5°: 300 mm

w wheel camber 3°: 350 mm

Folding height: 650 / 1000 / --- mm

Permitted inclination/slopes

max. obstacle height (depending on the setting of the leg support height):..... 0 to 100 mm
min. turning radius diameter:..... 1200 mm
Max. permissible rising gradient: 4.5° (8 %)
Max. permissible falling gradient: 4.5° (8 %)
Max. permissible transverse gradient: 4.5° (8 %)
static tilting safety in all directions: 6° (10 %)

Climatic data:

Ambient temperature: -25 °C to +50 °C
Storage temperature: -40 °C to +65 °C



Weights

Permissible total weight: max. 160 kg
max. user weight (including additional load): max. 150 kg
Max. additional loading: 10 kg
Empty weight: min. 10.5 / max. 18 kg

Heaviest single component:

Drive wheel: 2.2 kg
Transport weight: min. 6.5 kg
(without arm supports, cushion, drive wheels, leg supports)

Meaning of the labels on the wheelchair

	<p>Attention! Read the operating manuals and other provided documentation.</p>
 <small>Ident. No. 9 440 238</small>	<p>Anchor point for transport safety of the wheelchair.</p>

Meaning of the symbols on the type plate



Manufacturer



Order number



Serial number




Production date (Year – Calendar week)



Permitted user weight



Permitted overall weight



Permitted axle weights



Permitted rising gradient



Permitted falling gradient

max. ... km/h Permitted maximum speed



The product is approved as a seat within a motor vehicle



The product is **not** approved as a seat within a motor vehicle.

INSPECTION CERTIFICATE

Vehicle data:

Model:

Delivery note no.:

Serial number (SN):

Recommended safety inspection 1st year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 2nd year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 3rd year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 4th year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 5th year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

WARRANTY / GUARANTEE

We accept legal liability for this product within the scope of our general terms and conditions and warranty and the guarantee according to our described quality service. For warranty and guarantee demands please contact your specialist dealer with following Warranty/Guarantee section and the there included information on model description, delivery note number with delivery date and serial number (SN).

The serial number (SN) can be read off of the type plate.

Pre-condition for the acceptance of liability in any case is the intended use of the product, the use of original spare parts by authorised dealers as well as maintenance and inspections in regular intervals.

Guaranty is not granted for surface damages, tyres of the wheels, damages due to loosened screws or nuts as well as worn out attachment holes due to frequent assembly work.

Interferences through radiation sources such as mobile phones with high transmission power, HiFi-equipment and other extreme interference radiators outside of norm specifications cannot be declared as warranty or guarantee claims.

Attention:

- ! Failure to observe the instructions in the operating manual, improperly carried out maintenance work and, especially,

technical changes and additions (add-ons) carried out without our prior consent will lead to a general loss of guarantee and product liability.

Note:

This operating manual as a part of the product is to be handed out in case of a change of owner.

We reserve the right to make technical improvements.



The product conforms with the EC Directive 93/42/EEC (MDD) for medical products.

Warranty / Guarantee section

Please fill out! Copy if necessary and send the copy to the specialist dealer.

Warranty / Guarantee

Model designation:

Delivery note no.:

SN (view type plate)

Date of delivery:

Stamp of the specialist dealer:

Inspection certificate for transfer

Vehicle data:

Serial number (SN):

Model:

Delivery note no.:

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Your specialist dealer

MEYRA GmbH



Meyra-Ring 2
D-32689 Kalletal-Kalldorf



Tel +49 5733 922 - 311
Fax +49 5733 922 - 9311



info@meyra.de

www.meyra.de
