

# Mio Move All prices valid from February 12th, 2018





□ anodized Gold Orange (04)□ anodized Violet (05)

Nº 1570..

German Aid Index Number. 18.99.02.1020

Light wheelchair for children and adolescents with mechanical seat tilting from -5° to +35° in the physiological turning point, optional in 2 frame forms (Retro and Taurus) and 2 differently strong V-formed abductions, back adjustable in angle of  $78^{\circ}$ -126° and folds down to the front, firm wheel base, back height grows + 5 cm (optional: seat width + 2 cm , seat depth + 3 cm), stabalising rod, low friction wheels, cable brake built in to the side guard, nominal load max. 50kg.



	balising rod, low friction wheels, cable brake built in to the side guard, nominal load max. 50kg.			
	☐ Quotation ☐ Order		,	
User/ID number	Delivery address/comp	oany stamp	Billing address	
(Please fill in with block letters)			same as delivery address	
First and surname	Contact person			
Health insurance company	Date and signature			
Frame versions		-		Retro frame
Retro frame N°59730		Taurus frame N°5973		
abduction LLL 15-34 cm casters on one line with rear wheels	□ Slight abduction LLL 15-31 cm  casters narrower than rear wheels	Strong abduction LLL 15-34 cm casters on one line with rear wheels	☐ Slight abduction LLL 15-31 cm  casters narrower than rear wheels	Taurus frame
Frame- and user m	easures			cf. example on p.6
Seat width in cm (10.)  20 22 24 26 28 30 32 34  we measure: inner distance side guard to side guard	18   20   22   24   26   28   30   32   34   36   we measure: front edge sear to front edge back rod	Back height in cm (12)  25 30 35 40  45  (Back height growable by +5 cm we measure: top edge seat board edge back rod  Not selectable when Tarta back sys	d to top	Real lower leg length of the user (LLL):    cm     Cushion thickness (CT)    cm
Colours	-		-	

## Moderate metallic effect (except for black and Ultramarine) Finish alternatives Special colours (without additional charge, except for Nightblue) ☐ Shiny ☐ Zinc Yellow (02) ☐ Dark Green (15) Strong metallic effect: Dormant colours: ☐ Pure Orange (21) ☐ Crimson (04) ☐ Silver Metallic (36) $\square$ Apple Green Dormant (18) ☐ Pastel Blue (20) ☐ Dark Blue (16) ☐ Anthracite Metallic (22) ☐ Fire Red Dormant (45) □ Matt ☐ Brilliant Blue (08) ☐ Violet Dormant (50) Not possible with special ☐ Cream White (17) not metallic: Coarse structured surface: Dormant col. with additional charge colours! ☐ Girly Pink (19) □ Black (31) ☐ Coarse Black (27) ☐ Nightblue Dorm. (47) ☐ Ultramarine Blue (37) ☐ Chrome Yellow (06) Nº 1561 Design package color anodization (Colour examples for anodized colours can be found on our website in the area equipment) ☐ anodized Blue (01) The following parts are anodized: Restrictions: ☐ anodized Red (02) Caster forks, supporting plate for side guards, hand rims The hand rims are only available in the design package ☐ anodized Green (03) anodized colours by choosing aluminum hand rims. (not

Frame options						
Charles of the Control of the Contro		• •	Step tube			
☐ Modification set for growable seat width + 2 cm and seat depth up to 3 cm	☐ by 7,5 cm to the back, and adjustment of axis, side guards and brake  Attention: Return to standard configuration only possible by new side guards!	swing away installation adjustable	☐ Step tube			
N°597799	N°760177		N°7408			


hand rims)

ting.


posssible with stainless steel, Maxgrepp, PVC, double or nob

The design package ist only possible in the listed set-

seat width - SW seat depth - SD back height - BH customized product discount - CPD







You can find the colours under www.sorgrollstuhltechnik.de.





Seat system	Back system				
Alu seat plate	Back angle	Model 1	Model 2	Model 3	
Siological turning point tilta- ble from -5°to +35°	Adjustable in angle by faster from 81°, 90°, 102°, 114° oder 126°, folds down forward, incl. stabilising rod	□ 15 mm (015)	BH+10 cm; curve depth:  15 mm (115) 30 mm (130) 50 mm (150) BH+20 cm; curve depth:	BH+5 cm; curve depth:  15 mm (315)  30 mm (330)  50 mm (350)  BH+10cm; curve depth:	
Blocking set		only from seat width 30cm	☐ 15 mm (215) ☐ 30 mm (230)	☐ 15 mm (415) ☐ 30 mm (430)	
Mounting set to limit tilting; towards the front to max.:  ☐ 0° tilting  towards the back to max: ☐ +30°, ☐ +25°, ☐ +20°, ☐ +15° tilting			☐ 50 mm (250)  Attention: curve depth 50mm only from seat width 30cm	☐ 50 mm (450) BH+15 cm; curve depth: ☐ 15 mm (515) ☐ 30 mm (530) ☐ 50 mm (550) Attention: curve depth 50mm only from seat width 30cm	
N°192810		N°693	N°693	N°693	

Back extension		Back cover	Tarta® back system	Fitting set
for curved back	Restrictions	Adjustable	Adaptation	for solid back
□ extension up to 17.5 cm, retractable, with detachable head cushion  from SW 22 cm, weight reduction of approx. 1000g compared to a head rest, not possible with head rest (N° 6860)	-Model 1: from BH 35 cm -Model 3: from BH 30 cm only with version BH+5 cm -Fitting set N° 692598: from SW 30 cm	incl. back cushion (thickness 2 cm) Breathable Betex cover □ black (20) □ blue (21) □ red (22) □ grey (23) □ green (24) □ yellow (25) □ orange (27)	Adaptation holder for Tarta back system incl. Tarta back tubes Mounting possible from SW26 Order information for Tarta back system and instructions at www.tartadesign.it/en/distributori/	quired for mounting a solid back that is manufactured by the specialist trade  Curve depth 50mm only
N°686077			N° 611705	N°692598

Seat cushions (All covers with zippers and changeable)					
	Standard cushion	Anatomical cushion	Cover design	Incontinence cover	
tick if no cushion is requi-	comfortable seating 3 cm (03)	Anatomically formed, slight abduction, changeable cover, thickness of cushion measured at slightest point: front 5 cm, back 3 cm. <b>Customized product</b>	☐ black (20) ☐ blue (21) ☐ red (22) ☐ grey (23) ☐ green (24) ☐ yellow(25)	For <b>standard cushion</b> as chosen on the left	
	N° 80	N° 891002	N° 80	N° 86	
Additional covers for cushions as chosen on the right  ☐ 1 cover (01)  ☐ 2 covers (02)	Pressure relief cushion for medium-soft seating  3 cm (03) 5 cm (05)	Additional cover for anatomic al cushion as configurated  Customized product	Seat cushion as chosen, without cover	☐ For <b>anatomical cushion</b> as chosen on the left	
N° 82	N° 81	N° 891003	N° 998001	N° 891004	


Back cushions (All covers with zippers and changeable)						
	Standard cushion	Anatomical back cushion	Cover design			
☐ No back cushion (please tick if no cushion is required)		depr. depth:   3 cm (04)  5 cm (06)  Only in association with back depression.	Breathable cover			
		N° 8911	N° 8517/8911			
Additional covers for cushions as chosen on the right ☐ 1 cover (01) ☐ 2 covers (02)		☐ Additional cover for anatomic al back cushion as configurated Customized product	☐ Back cushion as chosen, without cover			
N° 85	N° 8517	N° 891110	N° 998501			

seat width - SW seat depth - SD back height - BH customized product discount - CPD









# Push handles Standard one hand push handles push bail height adjustable additionaly can be folded and height adjustable with synthetic hinge and height adjustable with metal hinge













integrated in activity side

guards

Leg support with foot rest (incl. calf belt, adjustable in height-, depth- and angle)					
Model 1	Model 2	Calf rest	Heel cups		
	Foot rest divided, adjustable in distance and height, steplessly hinged upwards up to 180° From SW 22 cm and LLL 23 cm, max. LLL 35 cm	for leg support 401905,  Mounting calf rest  left light double-sided	Heel cups, self installation, black  ☐ Size 1 (01) 7 x 11.5 cm ☐ Size 2 (02) 9 x 12.5 cm ☐ Size 3 (03) 10 x 15.5 cm		
	N° 401905		N° 4020		
	□ Locking mechanism at 90° for foot plate Model 2 (N°401905)	Possible from LLL 25 cm, foot plates are no longer hingeable backwards. From SW 26 cm	Heel edge ☐ for foot plate/s		
N° 4009	N° 4009	N° 403106	N°401904		
Push handles					
Standard	Single hand push handle	Push bow made of plastic	or metal		
extendable, cranked:  15 cm (05)  25 cm (06)	☐ Foldable with steel joint, extendable by 20 cm only up to SW30 cm and nominal load max. 40kg	justable in angle made of plastic up to SW 30cm	standable by 20 cm, adjustable in angle made of metal, from SW 32cm		
		at SW20/22 and BH25/30 only li- mited height adjustment possible	at BH25/30 only limited height adjustment possible		
N° 65	N°660202	N° 6534	N° 6532		
Brakes					
Bowden cable brake					
☐ Bowden cable brake,	1	Trum brake  Drum brake for standard-,profile- or lightweight			

Wheels						
Standard wheel	Profile wheel	Lightweight wheel	Lightweight wheel			
☐ All-round wheel	reduces the risk of finger in-		with drum brake hub, very lightweight, incl. aluminium handrim			
Additional weight of approx. 900g/wheel compared to lightweight wheels	Additional weight of approx. 800g/wheel compared to lightweight wheels	"	Weight reduction of up to 1,1 kg/wheel compared to standard drum brake wheels.			

wheels (N° 21..30/31/82), incl. brake lever for the accompanying person

If Bowden cable brake for the user is additionally required please tick also


Only at a wheel camber of 11° (N°201516)

N° 220036


Tyres						
Standard	Standard	Design	Performance			
Pneumatic: ☐ Schwalbe Right Run 1" (grey/black) (13)	Solid tyres:  Standard tyres  (grey) (11)	☐ Schwalbe Right Run 1" (blue/black) (14) ☐ Schwalbe Right Run 1" (red/black) (15) <b>Only</b>	Pneumatic:  ☐ Schwalbe Marathon Plus, black, pneumatic flat free tyres 1"  Only for 22" and 24" wheels			
		N° 2024	N° 2012			



seat width - SW seat depth - SD back height - BH customized product discount - CPD







# Seat heights All measures WITHOUT cushion!

All measures ±5 %

Wheel combination	Seat height measured without cushion from seat edge to floor with horizontal seat		
	Info: Mounting the wheels in the wrong borehole leads to a tilting of the wheelchair		
20" x 4" (=B2)	37,5 cm		
20" x 5" (=B1)	37,5 cm		
22" x 4" (=B2)	39,5 cm		
22" x 5" (=B1)	39,5 cm		
24" x 4" (=B2)	42 cm		
24" x 5" (=B1)	42 cm		

Caster wheels						
Standard PU	Standard solid rubber	Standard pneumatic	Design PU	Design solid rubber		







Wheel size and wh	eel camber	Handrims Mountin	g	
Wheel size	Wheel camber	at standard wheel	at lightweight wheel	Ergonomic handrims
Size in inches  20" recommended for SD 18-28  22" recommended for SD 22-32  24" recommended for SD 26-36	Camber 9°  Mounting set wheel camber 7°/11°, preassembled on 7° 11°	Mounting  ☐ tight ☐ wide  ☐ high ☐ standard  20" only with mounting "high"  Material:  ☐ Stainless steel (3001)	Mounting  ☐ tight ☐ wide ☐ standard ☐ high  Material: Alu, silver anodized.  combination standard/high not possible with light-	☐ handrim with ergonomic handrim profile, black anodized with built-in rubber grip  Only for 24" wheels
	N° 201516	☐ Aluminium (3075)	weight wheel drum brake	N°342424
Double handrim	N 201010			IV 042424
with standard hub	with drum brake hub	Opposite side	Distance	
with standard hub and double handrim, only up to 9° camber Mounting □ left □ right	with drum brake hub and double handrim Mounting □ left □ right		Distance of double handrim  ☐ standard 18mm (11)  ☐ tight 9mm (12)	Distance tight: Not in com- bination with hand rim co- ver or coated hand rim
N° 2614	N° 2615	N° 261435		
Handrim options				
Silicone	PVC	Ergo rubber inlay	Maxgrepp handrims	
Detachable silicone cover, black (soft surface)  for standard wheel(20-24)  for profile wheel (02)  for lightweight wheel(05)	☐ Handrim coated with black PVC, (Tetra Grip) Not for profile wheels	☐ Rubber inlay, 13mm wide  Only for profile wheels	☐ Allows pushing with the palm of hand (soft surface)  Only for standard wheels  and lightweight wheels	
N° 32	N° 33	13 mm	1	
for double handrims	for double handrims  Handrim cover for double handrim, both handrims (inside and outside), coated with black PVC  N° 262402  Handrim cover for opposing side of the double handrim, coated with black PVC		33 mm	
N° 3205	N° 262401	N° 3201	N° 34	
Spoke cover	Design Decors (Net price	es are not discountable)		
Transparent	Please specify the desired mo	odel of the corresponding man	ufacturer:	
☐ transparent (00)	□ www.blueprince.de	□ www.speichen-schutz.de	□ www.wheelchairgadgets. com	☐ Foil for curved back plate with the corresponding design as chosen
	child's name on spoke cover	child's name on spoke cover	child's name on spoke cover	Only with curved back plate
Nº 25	N°2525	N°2526	N°2527	N° 2527
Caster wheels				
Standard PU	Standard solid rubber	Design PU	Design solid rubber	
Polyurethane:  4" (2812) 5" (2805)	Solid rubber: ☐ 4" (2804) ☐ 5" (2809)	Polyurethane, transparent, with LEDs:  4" (2820) 5" (2834)	Solid rubber, black, with alu rim  4" (2835) 5" (2836)	
Caster forks		Steering options		

☐ Alu fork, coated in

frame colour

N° 2829

☐ track fixation (straight)

N° **9102)** N° 2851

(Not in combination with shock protection

☐ Alu fork, black anodized

 $\hfill\square$  Locking mechanism for track fixation

(Not in combination with shock protection


N°9102)

N° 2850




seat width - SW seat depth - SD back height - BH customized product discount - CPD













## Frame options **Electric systems** For all electric aids we re-**Outdoor Front End** Splash guard Sun shade **Bumper** commend casters of 5<sup>e</sup> ☐ facilitates significantly ac-☐ Splash guard for Outdoor ☐ Bumper, made ☐ Variable adjustable sun higher. Please regard the tive and passive driving on Front End, of rubber to shade in black with meapossible mounting specifiprotect furniture, soft or rough surfaces, incl. made of black surements 36x36 cm (LxH) cations of the manufactujack up aid. plastic doors, etc. Can be Frame form "Taurus" not possible with mounted up to AAT Systeme www.aat-online.de a seat width Restraint system. of 38 cm. Alber Systeme N° 599516 N° 59951212 N°9102 N°910300 www.alber.de

Therapy tables				
Transparent	Design			
Central installation, plugged under seat plate, adjustable in height and distance  Size 0 (10) (39 cm width, recommended for SW 18-24)  Size 1 (11) (46 cm width, recommended for 4 SW 26-30)  Size 2 (12) (59 cm width, recommended for SW 32-36)  Restrictions in combination of therapy table and abduction wedge: Height adjustment of therapy table is limited and therapy table can only be used in combination with abduction wedge holders	Number of design for plexiglas therapy table according to wheelchairgadgets.com			
N° 9209	N° 920099			

**Positioning aids** Here you will find a small selection to mark. You will find more positioning aids, corresponding mounting possibilities and further equipment in the separate order form "Positioning aids"

Order form "positioning aids", available in download area.	2-point pelvic belt	4-point pelvic belt	Standard head rests	Lightweight head rest
T MAN TO THE PARTY OF THE PARTY	□ with plastic fastener □ with metal fastener The size is determined by us according to the chosen wheelchair measurements. N°780502	□ with plastic fastener     □ with metal fastener  The size is determined by us according to the chosen wheelchair measurements.  N°780504		Lightweight head rest, adjustable with two conical joints, connection and removable cover included. (cover with same Betex-color as back cushion)  Size 1 (20) Pad size 1.10x200mm
More products for positio- ning, e.g.:	Chest-shoulder belt	Protection cap	Model 5, (25) x-shaped small	☐ Size 2 (21) Pad size 130x260mm
belts,     fixations and corresponding connections,     abduction wedges,     side truss pads,     head rests and more On the order form refer to this new wheelchair order/enquiry	☐ Chest-shoulder belt without zipper	☐ protection cap for metal fastener  Makes independent opening for the patient difficult.	□ Model 6, (26) x-shaped large	Weight reduction approximately 800-1000g to standard head rest systems.
	N°7810	N°78050210	N°6860	N°686009

Transport brackets	
Transport brackets	
☐ Transport brackets (4 pieces) for transportation in a handicapped accessible passenger car	CRASH TESTED
N°58507710	

Lateral support							
Holder	Pad	Extensions					
Holder adjustable, size 1, incl. linkage, without pad  ☐ left side (01) ☐ right side (02)	Attention: Cushions for size 1 can only be mounted on truss pads size 1.	for truss pad size 2 Attention: Cushions for size 2 can only be mounted on truss pads size 2.					
N°685901  Holder adjustable, size 2, incl. linkage, without pad  □ left side (01) □ right side (02)	1 pcs. 2 pcs.	1 pcs. 2 pcs.	for size 2 (5,5 cm) (0203) 1 pos. 2 pos.				
N°685902	N°685901	N°685902	N°6859				



seat width - SW seat depth - SD back height - BH customized product discount - CPD

# **Terms and conditions**

 Preface and Scope
 (1) Our terms and conditions apply exclusively. Any deviating, opposing, or complementing terms and conditions on the part of the business partner do not form integral part of contract unless explicitly approved of by us. There are no exceptions to the consent requirement, which also applies in cases where we conduct deliveries without restrictions while aware of the business partner's terms and conditions

terms and conditions.

(2) Agreements made in the individual case with the purchaser (including side-agreements, supplements, and alterations) always take priority over these terms and conditions. A contract or our confirmation - both in written form - are authoritative for the content of such arrangements.

(3) Any legally relevant declarations and notifications the purchaser submits to us after conclusion of the content of such arrangements.

of the contract (e.g. setting of deadlines, notice of defects, declarations of rescission or abate-ment) require the written form to be valid. (4) Any information on the application of legal regulations serves the mere purpose of clarification.

All regulations are legally binding without any such clarifications unless explicitly altered or excluded in these terms and conditions.

All prices are ex works, excluding insurance and packaging. The prices are net prices; the value-added tax is stated and charged separately according to the legally binding rate at the time of invoicing. The prices valid on the day of delivery are binding unless explicitly arranged otherwise.

## 3. Conclusion of Contract

SORG's offers are invitatio ad offerendum (invitations to treat). Ordering goods counts as an obligatory offer of contract on the part of the purchaser. Unless the order suggests otherwise, we have the right to accept this contract offer within two weeks time after reception. The acceptance of an offer of contract is effected by the confirmation of the order.

# 4. Terms of Payment

(1) Provided that no other agreements exist, payment is to be made by paying in advance – 50% at placing of order, 50% pre-delivery of goods. The purchaser bears any costs arising from transactions or discounting.

(2) If a 30-day payment deadline is exceeded, a delay rate of 8% on top of the respective base rate of the European Central Bank falls due without the need for prior notices. Beyond that, the purchaser is charged 3,- € for any reminder sent. We reserve the right to make use of the commercial due rate according to § 353 HGB. A higher damage for delay may be enforced via a suitable verification.

All current price lists are generally valid until March 31 of the respective year, or until new price lists are published. The latest version of the price list will always be provided on our homepage. Offers older than 3 months at that point expire and must be replaced with new offers and prices.

6. Time of Delivery and Delay We state the terms of delivery when accepting the order. The delivery date is non-binding and has mere guideline value - especially with regard to customised products. The terms for delay in delivery are defined by legal regulations. In any case, a writing notice issued by the purchaser is necessary for them to take effect.

7. Force Majeure (Act of God) (1) SORG is not to be held responsible for the following delays of delivery. This also applies if the (1) SORG IS not to be need responsible for the following delays of delivery. This also applies it also applies a sec causes occur with respect to legal representatives, auxiliary persons, or deliverers of SORG: Circumstances of force majeure as well as other extraordinary events that are unpredictable and unavoidable for, and not caused by, SORG, and that only occur after ordering or reserving, and that SORG could not have known about. Moreover, this applies in case of strikes, fire, floods, labour disputes, business disruptions, altered magisterial approbation and legal situation, and magisterial orders which are not due to operational risks.

(2) In the cases of paragraph 1, SORG is eligible to delay delivery/performance for the time of hin-drance as well as an appropriate start-up time. If SORG thus causes unacceptable complications for the purchaser, both SORG and the purchaser are eligible to withdraw from the contract. Circumstances that only lead to a temporal and hence acceptable delay of delivery and performance are excluded from this right. The legal rights of the client remain unaffected.

We generally dispatch uncovered by insurance and at the recipient's risk. We transfer the risk of coincidental loss or deterioration of the goods as well as delay on to the client at the moment of passing the goods on to the forwarder, bearer, or any other person/institute responsible for delivery. SORG will make any effort to respect the client's wishes in terms of path and kind of delivery. Any costs caused by this (including agreed carriage free delivery) are to be covered by the client. In the case of carriage free delivery conducted by SORG, the recipient is required to put forward any carriage costs for a refund. If the client wishes the dispatch to be delayed, any risk of loss or deterioration of the goods is carried by the recipient from the first day of readiness for dispatch. The goods dispatched are only insured upon special demand and at the cost of the recipient.

9. Damage through DeliveryAngelieferte Gegenstände sind, auch wenn sie unwesentliche Mängel Upon delivery, the goods need to be checked for obvious damage in the presence of the bearer. Any damage is to be documented in a written form and co-signed by the bearer. This protocol is to be forwarded to us immediately. Delayed reclamations, or reclamations not signed by the bearer,

10. Warranted Features, Receipt, and Technical Alterations
Delivered objects, even with minor damage, are to be accepted by the client, even if §433.paragraph 1 BGB is not fulfilled. Partial deliveries are permissible. Features only count as warranted if explicitly designated as such. A warranty of condition can only be asserted explicitly and in a written form. SORG reserves the right for technical alterations, especially for solutions of equal or higher value.

(1) SORG warrants for any defects of delivered new goods for 12 months as long as the purchaser is not the user (§§ 474, 13 BGB). Should the purchaser be the user then the warranty period is two years for new goods, for used goods 12 months. Should the purchaser also be an entrepreneur, then the liability for used items is void, as long as there is no claim for compensation according

(2) The purchaser's claim for defects presuppose that he has attended his obligation to give notice of defects on the legal grounds for his reasons of inquiry (§§ 377, 381 HGB). If defects are found during the inspection or at a later point, we are to be notified in writing without delay. The notice counts as not delayed if issued within two weeks. Posting the notice on time suffices for the keeping of the term. Independently from this obligation to inspect and give notice, the purchase

is to report obvious defects (including wrong and short delivery) within two weeks upon delivery (again, posting the notice on time suffices). If the purchaser misses the proper inspection and/or notice of defects, we do not warrant for the reported defects. If the delivered good is defected, we may choose to conduct supplementary performance either by repairing the defect (subsequent improvement) or by delivering flawless goods (compensation delivery). Our right to refuse the chosen method of supplementary performance on the respective legal grounds remains

(3) Our liability for defects expires when alteration or maintenance works are done by the client or a third-person-party without our prior written permission. Consequences of excessive strain, negligent or incorrect handling, and violent damage are excluded from the warranty for defects.

(1) With respect to compensation - regardless of the legal cause - SORG is only liable in case of malice and act of gross negligence. In case of slight negligence, SORG is only liable for damage of life, body, health, and for damage caused by the negligence of an essential contractual obligation (an obligation whose fulfilment make the proper realisation of the contract possible in the first place and in whose fulfilment make the proper realisation of the contract possible in the first place and in whose fulfilment the contract partner regularly trusts and can trust). In this case, however, our liability does not exceed the covering of predictable, typically occurring damage. (2) A free right to termination on the part of the purchaser (acc. to § 651, 649 BGB in particular) is precluded. Otherwise, the usual prerequisites and legal consequences apply.

13. Costs of Supplementary Performance Costs of supplementary performance cannot be charged on SORG's account.

14. Voluntary Withdrawal of Goods SORG is not obliged to take back any delivered end products or parts thereof. If SORG voluntarily decides otherwise in individual cases (especially for therapeutically altered requirements or deaths), SORG refunds a maximum of 75% of the delivery costs for unused parts whose delivery does not date back more than three months. There is no legal obligation concerning withdrawal in this respect. Deliveries which date back more than three months as well as sewed parts and customised products are not returnable. A copy of the delivery note or the bill serves as proof of

15. Reservation of Ownership
Goods delivered by SORG remain the property of SORG until any obligations from the parties'
business contract on the part of the purchaser have been redeemed. This also applies if SORG
has adapted individual or all claims into a running account which has been balanced and
acknowledged. The purchaser has the right to dispose of the goods in the proper course of
business. The purchaser hereby yields all claims to SORG which arise from the resale against
the purchaser or third parties. The purchaser is authorised to withdraw these claims even upon
conveyance. SORG's right to withdraw the claims remains untouched. However, SORG is obliged
to withdraw the claims as land, as the unwhords complice with his duty, of sorganized. not to withdraw the claims as long as the purchaser complies with his duty of payment and any other duties. SORG has the right to demand that the purchaser makes known the yielded any other duties. SORG has the right to demand that the purchaser makes known the yielded claims and their debtors, that he provides all information necessary for conveyance, and that he hands in all associated documents. As long and insofar the reservation of ownership is intact, the purchaser must not transfer any goods for securing or bailment without SORG's permission. The purchaser may process or alter the goods for the benefit of SORG exclusively. If processed with other goods which do not belong to SORG, SORG obtains co-ownership of these new goods in relation to the invoice value to the cost price of the other processed goods. If a third party interferes with SORG's ownership, the purchaser is to notify SORG immediately. The purchaser must not make any agreements with other purchasers which may exclude or limit SORG's rights. In particular, the purchaser must not make any agreements which destroy or limit the advance assignments to SORG. If the estimated value of the existing guarantees exceeds the claims to be ensured by more than 20%, SORG is insofar obliged to release by request of the purchaser.

# 16 Convright Protection

In Copyright Protection (If SORG conveys any images, drawings, indications of weight or measurements, technical documentations, or other product descriptions and documents (including electronic ones), the respective indications only apply if explicitly stated that they do. SORG reserves ownership and copyright of all these documents. They may only be made available to third parties upon SORG's explicit permission.

17. Place of Fulfilment, Place of Jurisdiction, and Severability Clause
(1) The law of the Federal Republic of Germany applies for these Terms and Conditions and all legal relations between us and our business partners. All international and supranational (contractual) legal systems, the UN Convention on Contracts for the International Sale of Goods in particular, are hereby excluded. Prerequisites and effects of the reservation of ownership acc. to no. 15 are subject to the right of the respective storage location of the good if the choice of law made in favour of the German law is illegitimate or ineffective. The only place of jurisdiction (also on an international basis) for any disputes arising directly or indirectly from the contractual relationship is, according to German Commercial Code, our office in Oberhausen-Rheinhausen, Germany. However, we also have the right to file a suit concerning the general jurisdiction of the

(2) If any clause in these Terms of Conditions should turn out to be ineffective, this does not affect the effectiveness of all other clauses.

Technical alterations and misprints reserved. As of: January 2017 SORG Rollstuhltechnik GmbH + Co. KG Benzstr. 3-5, 68794 Oberhausen-Rheinhausen / Germany Phone +49 7254 / 9279-0 Fax +49 7254 / 9279-10 in fo@sorgrolls tuhl technik. de, www.sorgrolls tuhl technik. de