



Corbac™ Adjustable Contour Back Support



▲ First, dial in the exact posterior contour with the Corbac's adjustment straps.

Then, refine the fit even further with the thoraco-lumbar-pelvic (TLP) supports. TLP supports provide precise adjustment of posterior and lateral contours.



Lightweight

At two pounds, the Corbac™ weighs little more than standard back upholstery.

Precision fit

The Corbac is easily adjusted to help improve postural control and alignment. It does not decrease the wheelchair's depth or width. The sitter experiences fitting in, not sitting on, the wheelchair.



Reduce heat and moisture buildup

The spacer fabric overlay promotes air circulation and wicks moisture from the back surface.


Versatile

No screws, nuts, or tricky hardware, just a quick and easy mounting system that takes seconds to assemble onto the wheelchair.

A cinch to clean

The overlay detaches easily and can be machine washed.

Ride Designs® Corbac™ Adjustable Contour Back Support



The unique spacer fabric material is used in the cover design for improved comfort and breathability.



Ride Corbac Sizes

Fits Wheelchair Widths		Corbac Height		Part Number
12"	30cm	12 - 14"	30 - 36cm	CB-1200
14"	36cm	12 - 15"	30 - 38cm	CB-1400
15"	38cm	13 - 16"	33 - 41cm	CB-1500
16"	41cm	14 - 17"	36 - 43cm	CB-1600
16"	41cm	17 - 20"	43 - 51cm	CB-1600T
17"	43cm	14 - 17"	36 - 43cm	CB-1700
18"	46cm	14 - 17"	36 - 43cm	CB-1800
18"	46cm	17 - 20"	43 - 51cm	CB-1800T
19 - 20"	48 - 51cm	14 - 17"	36 - 43cm	CB-1920
19 - 20"	48 - 51cm	17 - 20"	43 - 51cm	CB-1920T

▶ Note: For people at higher risk for skin breakdown, and/or more challenging postural control needs, consider the **Ride Java® Back** or the **Ride Custom Back**, the ultimate in postural support and skin care.



Ride Designs®
a branch of Aspen Seating, LLC

phone: 303.781.1633
toll-free: 866.781.1633
www.ridedesigns.com

© 2017, Ride Designs. 090-007-A
Patents: www.ridedesigns.com/patents