



Unit G Sark Business Park  
 Purcellsinch Industrial Estate  
 Kilkenny  
 R95 HK81  
 Email: info@momentumhc.ie  
 Tel: 056 7712394

# Ride Forward™ Cushion Order Form

<b>Therapist Name &amp; Address</b>	<b>Delivery Address</b>	<b>Client Details</b>
_____	_____	_____
_____	_____	_____
_____	_____	_____
Phone _____	Phone _____	_____
Fax _____	Fax _____	Date _____

## Ride Forward Cushion (E2607, skin protection and positioning wheelchair seat cushion)

	Forward Cushion Size (width x depth)	Part Number Complete	Fits Wheelchair Widths	Fits Wheelchair Depths	Quantity	Total
P	10 x 10" 25 x 25cm	FCD-1010	9/10" 23/25cm	9/10" 23/25cm		
P	10 x 12" 25 x 30cm	FCD-1012	9/10" 23/25cm	11/12" 28/30cm		
P	12 x 12" 30 x 30cm	FCD-1212	11/12" 28/30cm	11/12" 28/30cm		
P	12 x 14" 30 x 36cm	FCD-1214	11/12" 28/30cm	13/14" 33/36cm		
P	12 x 16" 30 x 41cm	FCD-1216	11/12" 28/30cm	15/16" 38/41cm		
	14 x 14" 36 x 36cm	FCD-1414	13/14" 33/36cm	13/14" 33/36cm		
	14 x 16" 36 x 41cm	FCD-1416	13/14" 33/36cm	15/16" 38/41cm		
	16 x 16" 41 x 41cm	FCD-1616	15/16" 38/41cm	15/16" 38/41cm		
	16 x 18" 41 x 46cm	FCD-1618	15/16" 38/41cm	17/18" 43/46cm		
	16 x 20" 41 x 51cm	FCD-1620	15/16" 38/41cm	19/20" 48/51cm		
	18 x 16" 46 x 41cm	FCD-1816	17/18" 43/46cm	15/16" 38/41cm		
	18 x 18" 46 x 46cm	FCD-1818	17/18" 43/46cm	17/18" 43/46cm		
	18 x 20" 46 x 51cm	FCD-1820	17/18" 43/46cm	19/20" 48/51cm		
	20 x 16" 51 x 41cm	FCD-2016	19/20" 48/51cm	15/16" 38/41cm		
	20 x 18" 51 x 46cm	FCD-2018	19/20" 48/51cm	17/18" 43/46cm		
	20 x 20" 51 x 51cm	FCD-2020	19/20" 48/51cm	19/20" 48/51cm		

**Subtotal Ride Forward Cushions** \_\_\_\_\_

## Additional Forward Items

Item	Part Number	Retail Price	Quantity	Total
Outer Cover Only -- specify Cushion width x length --	FCE-_____			
Ride CAM® Wedge Kit* -- specify Cushion width --	FCD-WK _____			

**Subtotal Additional Forward Items** \_\_\_\_\_

**TOTAL from all categories** \_\_\_\_\_

\* NOTE: Each kit contains two (2) Ride CAM wedges for 14 - 20" (36 - 51cm) width cushions only. Patents: www.ridedesigns.com/patents © 2015, Ride Designs. 090-031-A

**P** = Contoured and sized for pediatric application.

### Return Policy

Ride Forward Cushions may be returned if they are in their original, sealed, Ride plastic bag. All returns are subject to a 10% restocking fee.

### Have you considered a back support?



### The Corbac™ Adjustable Contour Back Support:

- lightweight
- breathable
- easy adjustment
- no hardware
- folds with wheelchair

E2611, General use wheelchair back cushion

### Ride CAM® Wedge Kit



Use Ride CAM Wedges **symmetrically** (shown above) to promote more upright posture and decrease ischial pressure or **asymmetrically** to help level a flexible pelvic obliquity.

Each kit contains two (2) Ride CAM wedges for 14 - 20" (36 - 51cm) width cushions only.